

## BECKMANN SELFRETAINING RETRACTOR


blunt; 4X5 teeth

33.0cm (13.2")



**SKU:** 070-691-32

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV273
<b>KLS Martin</b>	15-777-32
<b>AUTHORS</b>	<a href="#">BECKMANN</a>
<b>Medicon</b>	24.32.06
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV273
<b>KLS Martin</b>	15-777-32

## AUTHORS

[BECKMANN](#)

## Medicon

24.32.06

## SPECIFICATIONS

[blunt](#), [toothed](#)