

## BENNETT BONE LEVER

60mm; blunt


**SKU:** 190-580-60

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK180
<b>KLS Martin</b>	23-561-24
<b>AUTHORS</b>	<a href="#">BENNETT</a>
<b>Medicon</b>	77.12.06



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>Aesculap</b>	FK180
<b>KLS Martin</b>	23-561-24
<b>AUTHORS</b>	<a href="#">BENNETT</a>
<b>Medicon</b>	77.12.06