

## BENNETT BONE LEVER

45mm; blunt




**SKU:** 190-580-45

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>AUTHORS</b>	<a href="#">BENNETT</a>
<b>Medicon</b>	77.12.04



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>AUTHORS</b>	<a href="#">BENNETT</a>
<b>Medicon</b>	77.12.04