


## BERLIN TUBING CLAMP

straight; serrated

14.5cm (5.8")



**SKU:** 032-006-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>Aesculap</b>	AN908
<b>KLS Martin</b>	13-940-14
<b>AUTHORS</b>	<a href="#">BERLIN</a>
<b>Medicon</b>	17.73.14
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>Aesculap</b>	AN908
<b>KLS Martin</b>	13-940-14

**AUTHORS**

[BERLIN](#)

**Medicon**

17.73.14

**SPECIFICATIONS**

[serrated](#), [straight](#)