


BERLIN TUBING CLAMP

straight; serrated

14.5cm (5.8")



SKU: 032-006-14

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Tubing Clamps |
| Aesculap | AN908 |
| KLS Martin | 13-940-14 |
| AUTHORS | BERLIN |
| Medicon | 17.73.14 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Tubing Clamps |
| Aesculap | AN908 |
| KLS Martin | 13-940-14 |

AUTHORS

[BERLIN](#)

Medicon

17.73.14

SPECIFICATIONS

[serrated](#), [straight](#)