

BERNHARD TOWEL FORCEPS

straight; sharp

16.5cm (6.6")



SKU: 030-044-16

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF426
KLS Martin	13-921-16
AUTHORS	BERNHARD
Medicon	17.54.16
SPECIFICATIONS	sharp , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF426
KLS Martin	13-921-16

AUTHORS

[BERNHARD](#)

Medicon

17.54.16

SPECIFICATIONS

[sharp, straight](#)