


BIEMER-MÜLLER APPROXIMATOR

jaw opening 5.0mm; length 7.0mm

closing force 30-40g



SKU: 026-540-07

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
Aesculap	FD613
AUTHORS	BIEMER , MUELLER
Medicon	56.76.09

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
Aesculap	FD613
AUTHORS	BIEMER , MUELLER
Medicon	56.76.09