


**BIEMER-MÜLLER APPROXIMATOR**

straight; jaw opening 5mm; length 9mm  
 closing force 30-40g



**SKU:** 026-540-09  
**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

**ZUSÄTZLICHE INFORMATIONEN**

|         |                                                  |
|---------|--------------------------------------------------|
| UNIT    | <a href="#">pcs.</a>                             |
| TYPE    | <a href="#">Clips, Approximators</a>             |
| AUTHORS | <a href="#">BIEMER</a> , <a href="#">MUELLER</a> |
| Medicon | 56.74.09                                         |

**ZUSÄTZLICHE INFORMATIONEN**

|         |                                                  |
|---------|--------------------------------------------------|
| UNIT    | <a href="#">pcs.</a>                             |
| TYPE    | <a href="#">Clips, Approximators</a>             |
| AUTHORS | <a href="#">BIEMER</a> , <a href="#">MUELLER</a> |
| Medicon | 56.74.09                                         |