


BIEMER-MÜLLER APPROXIMATOR

angled; jaw opening 5mm; length 7.5mm

closing force 30-40g



SKU: 026-541-08

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
Aesculap	FD614
AUTHORS	BIEMER , MUELLER
Medicon	56.75.09

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
Aesculap	FD614
AUTHORS	BIEMER , MUELLER
Medicon	56.75.09