

BIEMER VESSEL CLIP

jaw opening 5,0mm; length 9.0mm

closing force 20-25g



SKU: 026-506-03

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
Aesculap	FD562
AUTHORS	BIEMER
Medicon	56.87.18



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
Aesculap	FD562
AUTHORS	BIEMER
Medicon	56.87.18