

## BILIARY STENT; DOUBLE-PIGTAIL

polyethylene with distal gold ring

8,5fr./7cm(Flaps); length: 10cm; single use



**SKU:** PK-P857DP

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)