

## BILIARY STENT; DOUBLE-PIGTAIL

polyethylene with distal gold ring

8,5fr./9cm(Flaps); length: 12cm; single use



**SKU:** PK-P859DP

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)