

## BILIARY STENT; DOUBLE-PIGTAIL SINGLE USE

polyethylene with distal gold ring

8,5fr./14cm (Flaps); length: 17cm



**SKU:** PK-P8514DP

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)