

## BILIARY STENT; DOUBLE-PIGTAILBILIARY STL

polyethylene with distal gold ring

12fr./6cm (Flaps); length: 9cm; single use



**SKU:** PK-P126DP

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)