


BIPOLAR CABLE 3.0M

Berchthold; Martin

bipolar System; Umax: 0.8 KVp

ATTENTION: OUR BIPOLAR RESECTOSOCOPY SYSTEM CAN BE USED ONLY WITH HF-UNITS WHICH HAVE THE POSSIBILITY OF BIPOLAR CUT IN SALINE SOLUTION !!!



SKU: 287-205-30BR

Categories: [Endoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables