


## BIPOLAR CABLE 3.0M FLAT-CABLE PLUG

international Models; 2-pins

bipolar System; Umax: 0.8 KVp

ATTENTION: OUR BIPOLAR RESECTOSOCOPY SYSTEM CAN BE USED ONLY WITH HF-UNITS WHICH HAVE THE POSSIBILITY OF BIPOLAR CUT IN SALINE SOLUTION !!!



**SKU:** 287-215-30BR

**Categories:** [Endoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>