


BIPOLAR CABLE 3.0M

Erbe ICC/ACC; Storz

bipolar System; Umax: 0.8 KVp

ATTENTION: OUR BIPOLAR RESECTOSOCOPY SYSTEM CAN BE USED ONLY WITH HF-UNITS WHICH HAVE THE POSSIBILITY OF BIPOLAR CUT IN SALINE SOLUTION !!!



SKU: 287-200-30BR

Categories: [Endoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--------------------------------|
| UNIT | pcs. |
| RANGE | HF-Cables |
| TYPE | Bipolar Cables |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--------------------------------|
| UNIT | pcs. |
| RANGE | HF-Cables |
| TYPE | Bipolar Cables |