

BIPOLAR CABLE; FLAT-CABLE PLUG

Valleylab; international Models; 5m

unit side Ø4,0mm; Umax: 0.8 KVp

Instrument-side



Unit-side



SKU: 287-215-50

Categories: [Endoscopy](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Cables](#)
TYPE [Bipolar Cables](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Cables](#)
TYPE [Bipolar Cables](#)