

## BIPOLAR CABLE; FLAT PLUG

Berchthold, Martin units; 5m

unit side round plug; Umax: 0.8 KVp

Instrument-side



Unit-side



**SKU:** 287-205-50

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>