


BIPOLAR CABLE; FLAT PLUG

US-Models; 3m; Umax: 0.8 KVp

unit side banana plugs 2x Ø 4.0mm

No CE mark! No labelling!



SKU: 287-210-30

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables