


## BIPOLAR CABLE; FLAT PLUG

US-Models; 3m; Umax: 0.8 KVp

unit side banana plugs 2x Ø 4.0mm

No CE mark! No labelling!



**SKU:** 287-210-30

**Categories:** [Endoscopy](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>