

## BIPOLAR CABLE; FLAT PLUG

Valleylab; international models; 3m

unit side Ø 4.0mm; Umax: 0.8 KVp

Instrument-side



Unit-side



**SKU:** 287-215-30

**Categories:** [Endoscopy](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [HF-Cables](#)  
**TYPE** [Bipolar Cables](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [HF-Cables](#)  
**TYPE** [Bipolar Cables](#)