

**BIPOLAR COAGULATION FORCEPS**

straight; blunt; Ø 1.0mm

15.0cm (6"); flat plug



**SKU:** 022-102-15  
**Categories:** [Open Surgery](#), [Forceps](#)

**ZUSÄTZLICHE INFORMATIONEN**

|                |                                                  |
|----------------|--------------------------------------------------|
| UNIT           | <a href="#">pcs.</a>                             |
| TYPE           | <a href="#">Bipolar Coagulation Forceps</a>      |
| Medicon        | 88.17.06                                         |
| SPECIFICATIONS | <a href="#">blunt</a> , <a href="#">straight</a> |



**ZUSÄTZLICHE INFORMATIONEN**

|                |                                                  |
|----------------|--------------------------------------------------|
| UNIT           | <a href="#">pcs.</a>                             |
| TYPE           | <a href="#">Bipolar Coagulation Forceps</a>      |
| Medicon        | 88.17.06                                         |
| SPECIFICATIONS | <a href="#">blunt</a> , <a href="#">straight</a> |