

BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 1.0mm

15.0cm (6"); flat plug



SKU: 022-103-15
 Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.08
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.08
SPECIFICATIONS	angled , blunt