

BIPOLAR COAGULATION FORCEPS

straight; blunt; Ø 2.0mm

15.0cm (6"); flat plug



SKU: 022-104-15
 Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.16
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.16
SPECIFICATIONS	blunt , straight