

BIPOLAR COAGULATION FORCEPS

straight; pointed; Ø 0.3mm

20.0cm (8"); flat plug



SKU: 022-100-20
 Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.20.05
SPECIFICATIONS	sharp , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.20.05
SPECIFICATIONS	sharp , straight