

BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 2.0mm

24.0cm (9.6"); flat plug

SKU: 022-105-24
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
 TYPE [Bipolar Coagulation Forceps](#)
 SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
 TYPE [Bipolar Coagulation Forceps](#)
 SPECIFICATIONS [angled](#), [blunt](#)