

BIPOLAR COAGULATION FORCEPS

angled on shanks; blunt; Ø 1.0mm

17.0cm (6.8"); flat plug

SKU: 022-107-10
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.26
SPECIFICATIONS	blunt , curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.26
SPECIFICATIONS	blunt , curved