

BIPOLAR COAGULATION FORCEPS

bajonett; straight; pointed; Ø 0.3mm

16.0cm (6.4"); flat plug



SKU: 022-110-16
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.67.05
SPECIFICATIONS	sharp , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.67.05
SPECIFICATIONS	sharp , straight