

BIPOLAR COAGULATION FORCEPS

bajonett; straight; pointed; Ø 0.3mm

18.0cm (7"); flat plug



SKU: 022-110-18

Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [sharp](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [sharp](#), [straight](#)