

BIPOLAR COAGULATION FORCEPS

bajonett; straight; pointed; Ø 0.3mm

22.0cm (8.8"); flat plug



SKU: 022-110-22
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.23.03
SPECIFICATIONS	sharp , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.23.03
SPECIFICATIONS	sharp , straight