

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

16.0cm (6.4"); flat plug



SKU: 022-114-16
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [straight](#)