

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

18.0cm (7"); flat plug



SKU: 022-114-18
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	blunt , straight