

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

22.0cm (8.8"); flat plug



SKU: 022-114-22
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	blunt , straight