

BIPOLAR COAGULATION FORCEPS

bajonett; straight; blunt; Ø 2.0mm

24.0cm (9.6"); flat plug



SKU: 022-114-24
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [blunt](#), [straight](#)