

BIPOLAR COAGULATION FORCEPS

bajonett; angled up.; pointed; Ø 0.3mm

22.0cm (8.8"); flat plug



SKU: 022-121-22
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [sharp](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [sharp](#)