

BIPOLAR COAGULATION FORCEPS

bajonett; angled up.; pointed; Ø 0.3mm

24.0cm (9.6"); flat plug



SKU: 022-121-24
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	angled , sharp



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	angled , sharp