

BIPOLAR COAGULATION FORCEPS

bajonett; angled up; blunt; Ø 2.0mm

20.0cm (8"); flat plug



SKU: 022-125-20
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)