

BIPOLAR COAGULATION FORCEPS

angled; pointed; Ø 0.3mm

15.0cm (6"); flat plug



SKU: 022-101-15
 Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.07
SPECIFICATIONS	angled , sharp



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.07
SPECIFICATIONS	angled , sharp