

BIPOLAR COAGULATION FORCEPS

bajonett; angled up; blunt; Ø 2.0mm

24.0cm (9.6"); flat plug



SKU: 022-125-24
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
SPECIFICATIONS	angled , blunt