

BIPOLAR COAGULATION FORCEPS

bajonett; angled down; blunt; Ø 1.0mm

20.0cm (8"); flat plug



SKU: 022-133-20
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.70.39
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.70.39
SPECIFICATIONS	angled , blunt