

BIPOLAR COAGULATION FORCEPS

bajonett; angled down; blunt; Ø 1.0mm

24.0cm (9.6"); flat plug



SKU: 022-133-24
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)