

BIPOLAR COAGULATION FORCEPS

bajonett; angled down; blunt; Ø 2.0mm

18.0cm (7"); flat plug



SKU: 022-135-18
Categories: [Open Surgery](#), [Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)