

## BLOOD PRESSURE MANOMETERS

calibrated; with velcro cuff



**SKU:** 100-570-00

**Categories:** [Diagnostics](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>                      |
| <b>TYPE</b>       | <a href="#">Blood Pressure Manometers</a> |
| <b>KLS Martin</b> | 17-972-00                                 |

### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>                      |
| <b>TYPE</b>       | <a href="#">Blood Pressure Manometers</a> |
| <b>KLS Martin</b> | 17-972-00                                 |