

## BLOUNT KNEE RETRACTOR

45mm; 2 teeth

26.0cm (10.4")

**SKU:** 190-576-26

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>KLS Martin</b>	23-565-26
<b>AUTHORS</b>	<a href="#">BLOUNT</a>
<b>Medicon</b>	77.12.52



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Levers</a>
<b>KLS Martin</b>	23-565-26
<b>AUTHORS</b>	<a href="#">BLOUNT</a>
<b>Medicon</b>	77.12.52