

## BÖHLER BONE RONGEUR

straight; double convertible

15.0cm (6")


**SKU:** 196-307-01

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Rongeurs</a>
<b>KLS Martin</b>	23-876-15
<b>AUTHORS</b>	<a href="#">BOEHLER</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bone Rongeurs</a>
<b>KLS Martin</b>	23-876-15
<b>AUTHORS</b>	<a href="#">BOEHLER</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>