

BONNEY DRESSING FORCEPS

straight; serrated

18.0cm (7.2"); width of jaw 3.2 +/-0.2

SKU: 020-120-18

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD078
AUTHORS	WAUGH
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD078
AUTHORS	WAUGH
SPECIFICATIONS	serrated , straight