

**BORCHARDT FLAT TREPHINE**

with eye; shaft 5mm  
 Ø 1.8mm



**SKU:** 190-284-18  
**Categories:** [Open Surgery](#), [Orthopaedics](#)

**ZUSÄTZLICHE INFORMATIONEN**

|            |                             |
|------------|-----------------------------|
| UNIT       | <a href="#">pcs.</a>        |
| TYPE       | <a href="#">Hand Drills</a> |
| KLS Martin | 22-194-18                   |
| AUTHORS    | <a href="#">BORCHARDT</a>   |

**ZUSÄTZLICHE INFORMATIONEN**

|            |                             |
|------------|-----------------------------|
| UNIT       | <a href="#">pcs.</a>        |
| TYPE       | <a href="#">Hand Drills</a> |
| KLS Martin | 22-194-18                   |
| AUTHORS    | <a href="#">BORCHARDT</a>   |