

BORCHARDT TWIST DRILL

short working end; shaft 5mm

Ø 1.0mm



SKU: 190-283-10
 Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-192-10
AUTHORS	BORCHARDT

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-192-10
AUTHORS	BORCHARDT