

BORCHARDT TWIST DRILL

long working end; shaft 5mm

Ø 5.0mm



SKU: 190-282-50
Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-190-50
AUTHORS	BORCHARDT

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-190-50
AUTHORS	BORCHARDT